

An *In Vitro* Recovery Study of Sodium Phenylbutyrate and Taurursodiol From 3 Types of Dosing Containers and Various Percutaneous Endoscopic Gastrostomy Feeding Tubes




Machelle Manuel, Shonna Scalfani, Jacky Lei, Shelley Capurso

Amylyx Pharmaceuticals, Inc., Cambridge, Massachusetts

BACKGROUND

- Sodium phenylbutyrate and taurursodiol (PB&TURSO) is approved to treat amyotrophic lateral sclerosis (ALS) in the United States¹ and approved with conditions in Canada for the treatment of ALS²
- PB&TURSO is considered an oral suspension, not an oral solution¹
 - According to the United States Pharmacopeia Expert Committee on Nomenclature and Labeling Nomenclature Guidelines, suspensions are defined as liquid preparations containing drug substance(s) and consist of solid particles dispersed throughout a liquid phase in which the particles are present in excess of the solubility³
- PB is a white or yellow powder that is freely soluble in water, and TURSO is an amphiphilic bile acid, a white microcrystalline powder that is sparingly soluble in water¹
- PB&TURSO combination is a white to yellow powder for oral suspension that consists of fine to large granules¹
- When prepared for administration after stirring vigorously in water (as prepared according to the recommended instructions), PB&TURSO may appear cloudy and/or contain residue^{1,4} (Figure 1)

FIGURE 1. PB&TURSO SUSPENSION



- The efficacy and safety of PB&TURSO combination were assessed in the CENTAUR clinical trial, a phase 2, multicenter trial in adults living with amyotrophic lateral sclerosis (ALS)⁵
- Trial participants (N=137) were randomized 2:1 to PB&TURSO (n=89) or matching placebo (n=48) administered once daily by mouth or via feeding tube for 3 weeks and then twice daily (1 packet in the morning and 1 packet in the evening), if well tolerated, for the 24-week randomized controlled phase⁵
 - Placebo was matched with PB&TURSO for appearance, taste and solubility
- For all participants, the study drug was stirred vigorously in 8 ounces (≈237 mL) of room temperature water; the mixture was either consumed orally or administered via feeding tube for participants with a gastrostomy or nasogastric tube within 1 hour of mixing⁶
- Twenty-eight CENTAUR participants were reported as using a feeding tube—16 participants (18%) receiving PB&TURSO and 12 participants (25%) receiving placebo^{4,5}
- The adverse event profile was similar between PB&TURSO and placebo participants using a feeding tube in the trial⁴

OBJECTIVE

- To evaluate PB&TURSO oral suspension administration via different feeding tube and dosing container combinations as part of an *in vitro* recovery study

METHODS

Feeding Tubes and Dosing Containers

- Five types of percutaneous endoscopic gastrostomy (PEG) feeding tubes, also known as G tubes (manufactured by Avanos Medical, Inc.) (Table 1), and three types of dosing containers were tested (Table 2)
 - Several combinations of feeding tubes and dosing containers were used to evaluate PB&TURSO recovery via chromatography and were visually inspected (Figure 2)

TABLE 1. FEEDING TUBES TESTED

Feeding Tube	1	2	3	4	5
Description	One-piece PEG long G-Tube with balloon	Low-profile PEG G-Tube with balloon with right-angle, 12-inch ENFit extension set	Low-profile PEG G-Tube with balloon with right-angle, 12-inch catheter tip extension set	Low-profile PEG G-Tube with retention bolster with right-angle, 12-inch ENFit feeding extension set	Low-profile PEG G-Tube with retention bolster with right-angle, 12-inch catheter tip extension set
Material of construction	Medical-grade silicone	Medical-grade silicone	Medical-grade silicone	Medical-grade silicone	Medical-grade silicone
Tube size, Fr	14	12	12	14	14
Stoma length, cm	24	4	4	3	3

TABLE 2. DOSING CONTAINERS TESTED

Dosing Container	A	B	C
Description	ENFit tip syringe used as dosing container comes with feeding tube kit	50-mL Luer-Lok™ tip syringe (used as catheter tip)	Gravity feeding bag including ENFit connector, PVC, ABS, PP material
Material of construction	Medical-grade silicone	Medical-grade silicone	Medical-grade silicone
Volume, mL	35	50	1000

ABS, acrylonitrile butadiene styrene; PP, polypropylene; PVC, DEHP-free polyvinyl chloride.

FIGURE 2. COMBINATION OF FEEDING TUBES AND DOSING CONTAINERS USED FOR SAMPLE ADMINISTRATION



In Vitro Recovery Assay

- An *in vitro* recovery study was conducted using 1 packet of PB&TURSO for the morning dose and 1 packet of PB&TURSO for the evening dose (each packet mixed with 8 ounces [≈237 mL] of room temperature water)
 - Each feeding tube setup was flushed with 20 mL of water before and after dose administration
- Two samples for each combination of feeding tube and dosing container were prepared and assayed via high-performance liquid chromatography to determine the percentage of PB&TURSO at feeding tube exit
 - Room temperature water was used as dosing media
 - All samples were mixed for up to 20 minutes or when reconstitution was achieved
 - Each test set was dosed 2 times to represent the morning dose and the evening dose
- The % recovery of PB&TURSO at the feeding tube exit was calculated as follows:

$$\% \text{ recovery} = \frac{\% \text{ LC}_{\text{feeding tube}}}{\% \text{ LC}_{\text{initial dispersion}}} \times 100$$

where:


- LC = concentration of PB or TURSO in the solution
- % LC feeding tube = % LC of PB or TURSO at feeding tube exit from assay sample solution
- % LC initial dispersion = % LC of PB or TURSO from control sample solution

- Acceptance criterion for *in vitro* recovery was defined as within 90% to 110% recovery of PB&TURSO at the feeding tube exit relative to the control

RESULTS

- All combinations of feeding tubes and dosing containers described met defined acceptance criteria (Figure 3)

FIGURE 3. PB&TURSO FEEDING TUBE RECOVERY^a



^aValues displayed on this figure reflect mean % recovery Dose, dosing container; tube, feeding tube. For more information about the feeding tubes and dosing containers tested, refer to table 1 and table 2, respectively.

- The feeding tubes and dosing containers used underwent visual examination, and no sedimentation, aggregation, adherence, build-up, or clogging were observed

CONCLUSION

In vitro data support that PB&TURSO can be used in a variety of PEG tube and dosing container combinations

Acknowledgments

The authors would like to thank Bharat Joshi, Maureen McLaughlin, and their team at ThermoFisher Scientific for their support with conducting the *in vitro* feeding tube studies.

This study is sponsored by Amylyx Pharmaceuticals, Inc.

Disclosures

MM, SS, JL, and SC are full-time employees of and may have stock ownership in Amylyx Pharmaceuticals, Inc.

References

- RELYVRIO. Prescribing Information. Amylyx Pharmaceuticals, Inc.; 2022.
- ALBRILOZA. Product Monograph. Amylyx Pharmaceuticals, Inc.; 2023.
- United States Pharmacopeia Expert Committee on Nomenclature and Labeling. Nomenclature Guidelines. 2020.
- Data on File. Amylyx Pharmaceuticals, Inc.
- Paganoni S, et al. *New Engl J Med.* 2020;383(10):919-930.
- Paganoni S, et al. Supplementary appendix. *N Engl J Med.* 2020;383(10):919-930. Accessed August 16, 2023. <https://www.nejm.org/doi/full/10.1056/nejmoa1916945>.